



HeloSpray Device Specifications

- Sizes: 3/8", 1/2", 3/4", 1", 1-1/2", 2"
- Materials are BPE, FDA, and USP Class VI compliant.
- Rotating Spray Patterns: 120 degree, 180 degree, 360 degree.
- The HeloSpray Device can be designed to meet the Customer's cleaning application.
- The HeloSpray Device can be designed to fit and clean a wide variety of tank sizes.
- Multiple rotating HeloSpray device installation is available to avoid shadowing from internal tank components.
- Reliable rotating solution coverage accomplished with minimal engineering time.
- Limited engineering required for developing rotating spray pattern to achieve Riboflavin testing.

Application

- Effective on standard dome and flat top tanks.
- Cascading coverage erodes residue from tank wetted surfaces.
- By replacing a standard sprayball with the rotating action of the Fluid Line Helospray, weeks of engineering time is saved during the tank manufacturing stage.

